MSwab®

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Cat N.	製品	数量	使用方法*
6E011N	1ml MSwab® transport and preservation medium in 12x80mm screw cap tube 	300 pieces 6 boxes of 50 pieces	
6E012N	1ml MSwab® transport and preservation medium in 12x80mm screw cap tube + 1 regular FLOQSwabs®  	300 pieces 6 boxes of 50 pieces	Nasal, throat, vaginal, groin, armpit, rectal, wound and faeces
6E013N	1ml MSwab® transport and preservation medium in 12x80mm screw cap tube + 1 thin & flexible FLOQSwabs®  	300 pieces 6 boxes of 50 pieces	Nasopharyngeal
6E076N	3 ml MSwab® transport and preservation medium in 12x80 screw cap tube 	300 pieces 6 boxes of 50 pieces	
6U019N	2ml MSwab® transport and preservation medium in 12x80mm screw cap tube 	300 pieces 6 boxes of 50 pieces	
6E092N01	3ml MSwab® transport and preservation medium in 12x80mm screw cap tube + 1 thin & flexible FLOQSwabs®  	300 pieces 6 boxes of 50 pieces	Nasopharyngeal

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- 2E013S50 - Tube with glass beads
- PFT1Wg13R100 - Closure white caps

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1. Schnee SV et al (2017) Performance of the Alere i RSV assay for point-of-care detection of respiratory syncytial virus in children. *BMC Infect Dis.* 17(1):767
2. Mahony J et al (2013) Development of a sensitive loop-mediated isothermal amplification assay that provides specimen-to-result diagnosis of respiratory syncytial virus infection in 30 minutes. *J Clin Microbiol.* 51(8):2696-701
3. Mahony J et al (2013) Multiplex loop-mediated isothermal amplification (M-LAMP) assay for the detection of influenza A/H1, A/H3 and influenza B can provide a specimen-to-result diagnosis in 40 min with single genome copy sensitivity. *J Clin Virol.* 58(1):127-31
4. Peters RM et al (2017). Evaluation of Alere i RSV for Rapid Detection of Respiratory Syncytial Virus in Children Hospitalized with Acute Respiratory Tract Infection. *Journal of clinical microbiology*, 55(4), 1032–1036.
5. Badman SG et al (2021) A comparison of ThinPrep against four non-volatile transport media for HPV testing at or near the point of care. *Pathology.* 53(2):264-266
6. Peterson LR et al (2017) Performance of the cobas MRSA/SA Test for Simultaneous Detection of Methicillin-Susceptible and Methicillin-Resistant *Staphylococcus aureus* From Nasal Swabs. *Am J Clin Pathol.* 148(2):119-127
7. von Allmen N et al (2019) Liquid and Dry Swabs for Culture- and PCR-Based Detection of Colonization with Methicillin-Resistant *Staphylococcus aureus* during Admission Screening. *Eur J Microbiol Immunol (Bp).* 9(4):131-137
8. Heers T ET AL (2017) Loop-mediated isothermal amplification (LAMP) assay-A rapid detection tool for identifying red fox (*Vulpes vulpes*) DNA in the carcasses of harbour porpoises (*Phocoena phocoena*). *PLoS One* 12(9):e0184349

This document may contain product information otherwise not accessible or valid in your country. Please be aware that Copan Italia S.p.A. does take any responsibility for accessing such information which may not comply with any valid legal process, regulation, registration or usage in the country of your origin. Product clearance and availability restrictions may apply in some Countries. Please refer to Copan website (www.copangroup.com) to view and/or download the most recent version of the brochure. This document is mainly intended for marketing purposes, always consult product insert for complete information. The use of these products in association with diagnostic kits or instrumentation should be previously validated by the user. ©2021 Copan Italia. All rights reserved. The trademarks mentioned herein are property of Copan Italia S.p.A.
Code: JMKPF05R01.JA

Copan Italia s.p.a.
Via Francesco Perotti 10,
25125 Brescia, Italy

t | f +030 2687211
@ | info@copangroup.com
www.copangroup.com